

Multienzyme Isothermal Rapid Amplification

Basic kit (DNA)

Guide Manual

(Part number: WLB8201KIT)

Overview of the MIRA technology

Multienzyme isothermal rapid amplification(MIRA) is based on a combination of polymerases and DNA recombination/repair proteins: Recombinases bind to single-stranded nucleic acid and stimulate the resulting protein-DNA complex to search for homologous sequences in duplex DNA with Single Strand Binding protein. Once homology is located, a strand-switching reaction is performed and the oligonucleotide is paired to its complement permitting a polymerase to begin synthesis from the 3' end. The amplification process is very rapid when optimised and can reach detectable levels of product in a few minutes.

Primer design considerations

The recommended primer length is 30 to 35 nucleotides . It is best to avoid unusual sequence elements within the primer, such as long tracks of one particular nucleotide or a large number of small repeats. oligonucleotides that contain sequence elements that promote secondary structures and primer-primer interactions or hairpins should be discarded.The recommended amplicon length is not exceed 500 bp, and ideally is between 150 -300 bp .

Materials provided

| Materials provided | Content |
|-------------------------------|-----------------|
| A buffer | 1 tube × 1.6 mL |
| B buffer | 1 tube × 150 μL |
| Positive control DNA template | 1 tube × 30 μL |
| Positive control primer Mix | 1 tube × 60 μL |
| Reaction unit | 48 |
| Guide manual | 1 |

Storage conditions

- Transport temperature: ≤ 20 °C constant temperature environment.
- Store at -20 °C without light, no repeated freezing and thawing.
- Full activity is guaranteed for 14 months.

Set up

1. Add 29.4 μL A buffer (supplied) to a MIRA basic reaction.
2. Add 2 μL forward primer and 2 μL reverse primer (Primer concentration: 10 μM).
3. Add template and water (template and water total volume is 14.1 μL).

4. Add 2.5 μ l B buffer (supplied) and mix well, then spin briefly to start reaction.
5. Incubate at 37-39°C for 30 minutes.
6. After 30 minutes, clean amplicons before running on 1.5% agarose gels.

Reaction mix

| | |
|-----------------------------|--------------|
| A buffer | 29.4 μ L |
| Forward primer (10 μ M) | 2 μ L |
| Reverse primer (10 μ M) | 2 μ L |
| DNA template and water | 14.1 μ L |
| B buffer | 2.5 μ L |
| Total volume | 50 μ L |




Fig.1 M:DL2000 Marker, P: positive control, N: negative control.

Notes

1. Reactions start as soon as B buffer is added.
2. Set up a negative control for each experiment.
3. For research and development use only.

Ordering information and technical support

Amp-Future (Changzhou) Biotech Co., Ltd.
 4th Floor, No.9 Building, China Israel Changzhou Innovation Park (CICP),
 Changzhou City
 Jiangsu Province
 China