

Multienzyme Isothermal Rapid Amplification

RT- Basic kit (RNA)-II

Guide Manual

(Part number: WLRB8207KIT)

Overview of the MIRA technology

Multienzyme isothermal rapid amplification(MIRA) is based on a combination of polymerases and DNA recombination/repair proteins: Reverse transcriptase generate DNA from RNA template. Recombinases bind to single-stranded nucleic acid and stimulate the resulting protein-DNA complex to search for homologous sequences in duplex DNA with Single Strand Binding protein. Once homology is located, a strand-switching reaction is performed and the oligonucleotide is paired to its complement permitting a polymerase to begin synthesis from the 3' end. The amplification process is very rapid when optimised and can reach detectable levels of product in a few minutes.

Primer design considerations

The recommended primer length is 30 to 35 nucleotides . It is best to avoid unusual sequence elements within the primer, such as long tracks of one particular nucleotide or a large number of small repeats. oligonucleotides that contain sequence elements that promote secondary structures and primer-primer interactions or hairpins should be discarded.The recommended amplicon length is not exceed 500 bp, and ideally is between 150 -300 bp .

Materials provided

Materials provided	Content
A buffer	1 tube × 1.6 mL
B buffer	1 tube × 150 µL
Reaction unit	48
Guide manual	1

Storage conditions

- Transport temperature: ≤ 20 °C constant temperature environment.
- Store at -20 °C without light, no repeated freezing and thawing.
- Full activity is guaranteed for 14 months.

Set up

1. Add 29.4 µL A buffer (supplied) to a MIRA basic reaction.
2. Add 2 µL forward primer and 2 µL reverse primer (Primer concentration: 10 µM).
3. Add template and water (template and water total volume is 14.1 µL).
4. Add 2.5 µl B buffer (supplied) and mix well, then spin briefly to start reaction.

5. Incubate at 42°C for 30 minutes.
6. Add Tris saturated phenol/chloroform/isoamyl alcohol (25:24:1), and mix 1:1 to reaction solution, centrifuge at 12000 rpm for 5 min, take 5 µL of supernatant for 1.5% -2% agarose gels.

Reaction mix

A buffer	29.4 µL
Forward primer (10 µM)	2 µL
Reverse primer (10 µM)	2 µL
RNA template and water	14.1 µL
B buffer	2.5 µL
Total volume	50 µL

Notes

1. Reactions start as soon as B buffer is added.
2. Set up a negative control for each experiment.
3. For research and development use only.

Ordering information and technical support

Amp-Future (Changzhou) Biotech Co., Ltd.

4th Floor, No.9 Building, China Israel Changzhou Innovation Park (CICP),

Changzhou City

Jiangsu Province

China